



Restorative Care of America
INCORPORATED

Phone: (727) 573-1595 or (800) 627-1595
Fax: (727) 573-1886 or (800) 545-7938
12221 33rd Street N, St Petersburg, FL 33716
www.rcai.com

03/14 Literature may not reflect current manufacturer specifications.



LSO and TLSO with Universal Sizing



The new **Lumbar Sacral Support** and **Thoracic Lumbar Support** by RCAI feature “two sizes fits all”, cutting down on inventory and assuring proper patient fit. Each size offers two additional side panels, 3" and 6" to increase the circumference of the brace. The lightweight design of the LSO and TLSO offers a low profile approach in the stabilization of the lumbar sacral spine or thoracic spine following injury, post surgery or disease states. They feature an easy to apply design, including “back pack type” shoulder straps on the TLSO with patient friendly adjustments. Rigid, heat moldable Kydex® posterior/anterior panels and unique strapping design relieves load on intervertebral discs by lifting and holding the patient’s weight, improving posture and relieving pain. Bendable side panels with metal stays conform to the patient’s anatomy for comfort and support, and can be removed from the brace in the latter stages of the healing process.



The new Lumbar Sacral Support and Thoracic Lumbar Support feature heat moldable Kydex® anterior and posterior panels.

SIZING

	Measure Hip Circ.	LSO Model #	TLSO Model #
Size 1	33 – 48" (83.8 - 121.9 cm)	453LSO-S-S/L	463TLS-S-S/L
Size 2	48 – 62" (121.9 - 157.5 cm)	453LSO-S-XL/XXXL	463TLS-S-XL/XXXL

Limited Warranty Coverage: Any defect in material or workmanship solely limited to the Kydex® plastic frame or metal portion or hinge mechanisms for 1 (one) year from date of purchase.