



Phone: (727) 573-1595 or (800) 627-1595
 Fax: (727) 573-1886 or (800) 545-7938
 12221 33rd Street N, St Petersburg, FL 33716
 www.rcai.com

05/16 Literature may not reflect current manufacturer specifications.



Humeral Arm Brace System 94-HABS



Optional Supination/Pronation Attachment
 Specify Right or Left
Model # Right **Model # Left**
 74S/P-R 74S/P-L



Auxilliary Hinge
 Available for bilateral stability.
 Sold separately.
Model # 94-AUX-HNG

The new **Humeral Arm Brace System** by RCAI offers quadrilateral compression of the humerus and forearm providing optimal ROM and rotational control of the arm following post surgical or post injury indications. The hinge features flexion/extension ROM control setting of the elbow in 10° increments. An optional supination/pronation attachment is available for prescribed positioning. Manufactured of trimmable polyethylene, and lined with closed-cell foam, it is lightweight, durable, and cannot harbor moisture. Humeral and forearm stockinette included.

SIZING

Specify Right or Left

Size	Mid. Humeral Circ.	Model # Right	Model # Left
Small	7 1/2 – 9 1/2" (19.1 - 24.1 cm)	94-HABS-S-R	94-HABS-S-L
Medium	9 1/2 – 11 1/2" (24.1 - 29.2 cm)	94-HABS-M-R	94-HABS-M-L
Large	11 – 15" (27.9 - 38.1 cm)	94-HABS-L-R	94-HABS-L-L
XLarge	13 – 18" (33 - 45.7 cm)	94-HABS-XL-R	94-HABS-XL-L

Limited Warranty Coverage: Any defect in material or workmanship solely limited to the Kydex® plastic frame or metal portion or hinge mechanisms for 1 (one) year from date of purchase.