





# Introducing the award-winning Panthera CFII vertical shock foot

Continuous Single carbon spring with no interruption of bolts or laminated connections for comfortable, smooth and dynamic gait.

- Exceptional energy return
- Bumper-Carbon synergy provides dynamic assistance into midstance
- Continuous energy flow through single spring = no dead spots
- Lowest profile vertical-shock foot on the market
- Additional room for componentry
- Excellent ground compliance via freely suspended upper spring
- Individualized performance results via trimmable plug and play inserts
- Includes two footshells with integrated textile = no spectra socks required



*Winner of AVK design award for innovative elastomer/carbon integration.*

**Check out our CFII videos!**

[medi video showing athlete use outdoors](#)

[CFII feature video](#)

[CFII adjustment tips](#)