

# RUSH™ ROGUE

## HIGH IMPACT. EXTREME COMFORT.

RUSH™ ROGUE provides the most realistic and dynamic foot and ankle motion available. The Vertical Loading Pylon provides comfortable relief from high impact loading allowing the user to push the foot to the limits. The components are made of advanced fiberglass composite which is three times more flexible than most conventional prosthetic feet. The roll through characteristics of the foot provide exceptional energy return with no “dead spot” allowing users to push themselves to the limits on the court or in aggressive terrain.

- **Vertical Loading Pylon** for vertical shock relief and impact absorption
- **Smooth roll through** for high energy return and comfortable in any situation
- **Biomimetic ankle action** to simulate the human foot – no one will know you wear a prosthetic
- **Extreme reliability** for extended use in any high impact application a human can survive
- **Rotational torsion relief** to ease the impact to knee, hip and back

### PROFILE INFORMATION

<b>Weight Limit</b>	360lbs./150kg	<b>Heel Rise</b>	3/8"/10mm	<b>Size</b>	<b>Build Height*</b>
<b>Size</b>	22cm - 29cm	<b>†Suggested Codes</b>	L5987, L5984	<b>22 - 24cm</b>	6.5"
<b>Foot Shell</b>	Light or Dark	<b>Warranty</b>	36 Months	<b>25 - 27cm</b>	6.8"
<b>Product Weight</b>	865g/2.18lbs. (26cm)			<b>28 - 29cm</b>	7.3"

† The listing of this code should not be considered a guarantee for coverage or a payment. \* Build height includes footshell.

### CATEGORY CHART

Weight (lbs)	90 - 116	117 - 146	147 - 180	181 - 218	219 - 261	262 - 310	311 - 360
<b>Low Activity</b>	1	2	3	4	5	6	7
<b>Moderate Activity</b>	1	2	3	4	5	6	7
<b>High Activity</b>	2	3	4	5	6	7	8
<b>High Impact</b>	2-3	3-4	4-5	5-6	6-7	7-8	8-9

NOTE: Increasing category numbers indicate increased keel response.

SMOOTH · DYNAMIC · RELIABLE · VIRTUALLY INDESTRUCTIBLE