



Product Code

RUSH Kid = **KID**

Foot Shell = **FSK**

Spectra Sock = **SSK**

Activity Level

Moderate = **M**

High = **H**

High Impact = **P**

Age

Under 6 yrs = **01**

7 yrs to 9 yrs = **02**

10 yrs to 12 yrs = **03**

13 yrs and older = **04**

Foot Shell

Left Light = **L1**

Right Light = **R1**

Patient Weight

30 - 39 lbs = **03**

40 - 49 lbs = **04**

50 - 59 lbs = **05**

60 - 69 lbs = **06**

70 - 79 lbs = **07**

80 - 89 lbs = **08**

90 - 99 lbs = **09**

100 - 109 lbs = **10**

110 - 119 lbs = **11**

120 - 129 lbs = **12**

130 - 139 lbs = **13**

140 - 149 lbs = **14**

150 - 159 lbs = **15**

160 - 169 lbs = **16**

170 - 179 lbs = **17**

Foot Size

19 cm = **19**

20 cm = **20**

21 cm = **21**

Category Chart

Weight lbs	35 - 56	57 - 78	79 - 101	102 - 125
Moderate Activity	1P	2P	3P	4P
High Activity	1P	2P	3P	4P
High Impact	2P	3P	4P	5P

*If not selected, the stiffness will default based on weight and activity.

Activity Level Definitions:

Moderate - Jogging, Walking, Daily Activities, Golf, Skiing

High - Running, Baseball, Wakeboarding

High Impact - Sprinting, Basketball, Extreme Sports

* Customers ordering Foot Shell only will indicate: **Product Code, Foot Size and Foot Shell**

** Customers ordering Spectra Sock only will indicate: **Product Code and Foot Size**

Part Number Example: An **11 year old** amputee weighing **105 lbs** wears a **21 cm** foot and is considered a **high impact user**. The prosthetist would like to purchase a **RUSH Kid** foot with a **Right Light Foot shell**.



The Prosthetist would use the following part number nomenclature to place order:

