

mediven® for women vitality

compression your way



Features:

- Herringbone-ribbed design for a classic look at work or after work
- Luxury and style without sacrificing the medical accuracy you expect from medi
- Clima Fresh system inhibits bacteria growth and unpleasant odors
- Clima Comfort system combined with lanolin keeps you naturally dry and refreshed
- Available in all three compression classes: 15-20, 20-30, and 30-40 mmHg
- Available in standard and petite lengths

Available Options



- Closed Toe
- Standard Length
- Petite Length
- Chocolate
- Charcoal
- Ebony

Leg Circumference		I	II	III	IV	V	VI	VII
Calf	11"–14"	11¾"–15¼"	12¾"–16½"	13¾"–17¾"	14½"–19"	15½"–20¼"	16¼"–21½"	
	28–35.5cm	30–39cm	32.5–42cm	35–45cm	37–48.5cm	39.5–51.5cm	41.5–54.5cm	
Ankle	7"–7¾"	8"–8¾"	9"–9¾"	10"–10¾"	11"–11¾"	12"–13"	13¼"–14¼"	
	18–20cm	20.5–22.5cm	23–25cm	25.5–27.5cm	28–30cm	30.5–33cm	33.5–36cm	
Leg Length		Standard: 15"–19" (38–47.5cm)			Petite: 12"–15" (31–38cm)			

Compression that complements:
your style. your life. you.



* Not available in tan.

Introducing mediGrip®!

A revolutionary top band for added security.



Herringbone-ribbed pattern

mediGrip® is a completely seamless and integrated silicone top band for a calf-high stocking.

- No sewing seams - completely seamless
- Silicone thread is knit into the inside of the top band
- All the grip of a silicone top band without the irritations
- mediGrip® is skin-side only - no clinging clothing
- No compromise on efficacy - maintains gradient compression profile
- Standard feature on every pair of mediven® for women vitality

