

ENCORE II™

GEL LINER



- **Anti-elongation matrix**
- **Soft comfortable gel**
- **General purpose**

The Encore II is a soft general purpose liner designed for comfort and everyday use. The bottom third of the distal end is laminated with an anti-elongation matrix that is designed to reduce stretch and eliminate pistoning. Soft thermoplastic gel enriched with medical grade mineral oil ensures continuous comfort. The combination of advanced fabric design and gel results in a liner that provides reliable comfort and fit.

Thermoplastic Gel	Type	Color	Thickness	Size
Highly responsive	Locking Liner	Black	6mm/3mm	16cm 24cm
Motion control	Cushion Liner	Tan	9mm/3mm	18cm 28cm
Designed for comfort		Beige	6mm/6mm	20cm 34cm



(813) 606-4220
www.thermoplygel.com