



BREG LUMBAR BRACE FAMILY

INTRODUCING BREG'S NEW SPINE BRACES

Breg's goal is to provide support at every step of the orthopedic care journey to our customers and their patients. When it comes to spine products, that means making patient comfort and compliance top priorities. We've designed our new Pinnacle and Ascend™ lines with premium breathable fabrics and easy-to-adjust features, so patients get exactly the support they need, when they need it most. And since we manufacture these products in house, we are prioritizing patient experience in every detail.

All of our lumbar braces have patient-centric design features, including:

- Lightweight and breathable Airmesh® padding
- Universal sizing that makes them easy to fit
- 3:1 mechanical advantage lacer system
- Dual-pull lacer system with independent top and bottom adjustment
- Ergonomic pull handles

ASCEND LUMBAR ORTHOSIS

Our Ascend lumbar braces are designed as a universal option for stability and support of the lumbar spine.



Features of our Ascend lumbar orthosis include:

Lightweight and breathable
Airmesh padding



Dual-pull lacer system
with independent top
and bottom adjustment

Ergonomic pull handles

Ascend LO 627/642

SP10627-000



Ascend LSO 631/648

SP10631-000




+ Additional support
up to T9 vertebra

Size	1	2	3	4	5	6	EXT.
Waist Size	66.0–81.3cm 26" – 32"	78.7–94.0cm 31" – 37"	91.4–109.2cm 36" – 43"	106.7–124.5cm 42" – 49"	121.9–139.7cm 48" – 55"	137.2–152.4cm 54" – 60"	152.4cm+ 60" +

Ascend LSO 637/650

SP10637-000



+ Includes lateral pads
(provided with brace)

Size	1	2	3	4	5	6	EXT.
Waist Size	61.0–76.2cm 24" - 30"	73.7–88.9cm 29" - 35"	86.4–104.1cm 34" - 41"	101.6–119.4cm 40" - 47"	116.8–134.6cm 46" - 53"	132.8–147.3cm 52" - 58"	147.32cm+ 58" +

Ascend SI Belt 621

SP10621-0X0



Common examples of use Ascend SI Belt

- Sprain / Strain
- Facet Syndrome
- SI dysfunction
- Sacroiliitis
- Injections

Size	1	2	3
Waist Size	64 – 84cm 25" – 33"	84 – 124cm 33" – 49"	104 – 165cm 41" – 65"

ASCEND TLSO LUMBAR ORTHOSIS

Our Ascend TLSO brace is designed as a universal option for stability and support of the lumbar and thoracic spine.




Features of our Ascend TLSO orthosis include:



Ascend TLSO 456/457

SP10456-000



Common examples of use for TLSO braces

- Post-operative Support
- Fracture Management
- Decompressive Procedures
- Degenerative Disc Disease
- Kyphosis
- Scoliosis
- Kyphosis Secondary to Osteoporosis
- Spondylolisthesis
- Spinal Stenosis
- Spondylosis
- Bulging / Herniated Disc
- Severe DJD

Size	1	2	3	4	5	6	EXT.
Waist Size	66.0–81.3cm 26" – 32"	78.7–94.0cm 31" – 37"	91.4–109.2cm 36" – 43"	106.7–124.5cm 42" – 49"	121.9–139.7cm 48" – 55"	137.2–152.4cm 54" – 60"	152.4cm+ 60" +

