



BREG CERVICAL COLLAR FAMILY

INTRODUCING BREG'S NEW SPINE BRACES

Breg's goal is to provide support at every step of the orthopedic care journey to our customers and their patients. When it comes to spine products, that means making patient comfort and compliance top priorities. We've designed our new Pinnacle and Ascend™ lines with premium breathable fabrics and easy-to-adjust features, so patients get exactly the support they need, when they need it most. And since we manufacture these products in house, we are prioritizing patient experience in every detail.

Common Examples of Use for Cervical Collars

- Post-surgical Stabilization
- Cervical Disc Syndrome
- Cervical Stenosis
- Fracture Management
- Sprains / Strains
- Decompressive Procedures
- Cervicogenic Headaches
- Radiculopathy
- Trauma / Whiplash

PINNACLE CERVICAL COLLARS

Breg's new Pinnacle cervical collars have universal sizing for easy inventory management, plus patient-centric features to improve compliance.



Features of our Pinnacle cervical collars include:



The Pinnacle family includes three types of collars, so providers can choose the best product for their patients.

Pinnacle Cervical Collar 172

SP40172-000



Pinnacle Cervical Collar TX 174

SP40174-000



- + Anterior and posterior thoracic extensions provide additional support and control

Pinnacle Cervical Collar MP 180

SP40180-000



- + Anterior thoracic extension
- + Multi-Post back panel with pivoting occipital support plates that contour to the patient's head
- + Support plate positions can be adjusted to avoid discomfort on the occipital protuberance and nerves

ASCEND CERVICAL COLLARS

Ascend collars come in four adult sizes and five pediatric sizes to accommodate a wide range of patients and anatomical differences.



Features of our Ascend cervical collars include:

Mix and match front and back panels separately to accommodate unique anatomies

Specially shaped peripheral support tabs provide increased surface contact area with reduced flexing zone for enhanced comfort



Antimicrobial cotton-hemp blend liners to help prevent skin irritation

Ultra-breathable foam padding to reduce perspiration



An extra set of pads to allow for continuity of care during pad cleaning


Ascend Cervical Collar 172

SP30172-0X0



Ascend Pediatric Collar

SP30272-0X0



Size	PD1	PD2	PD3	PD4	PD5
Age	1-18 Months	9-24 Months	1-3 Years	2-5 Years	3-6 Years
Weight	5-13 kg 11-29 lbs	10-15 kg 22-33 lbs	11-16 kg 24-36 lbs	12-19 kg 26-42 lbs	12-25 kg 27-54 lbs
Length	53-84 cm 21-33 in	74-94 cm 29-37 in	84-102 cm 33-40 in	89-114 cm 35-45 in	94-122 cm 37-48 in

PINNACLE

Part Numbers	
Collars	
Pinnacle Cervical Collar 172	SP40172-000
Pinnacle Cervical Collar TX 174	SP40174-000
Pinnacle Cervical Collar MP 180	SP40180-000
Pad Kits	
Pad Kit, Pinnacle Cervical Collar 172	SP40272-000
Pad Kit, Pinnacle Cervical Collar TX 174	SP40274-000
Pad Kit, Pinnacle Cervical Collar MP 180	SP40280-000

ASCEND

Part Numbers	
Adult Collars	
Ascend Cervical Collar 172, Short	SP30172-020
Ascend Cervical Collar 172, Regular	SP30172-030
Ascend Cervical Collar 172, Tall	SP30172-040
Ascend Cervical Collar 172, X-Tall	SP30172-050
Pediatric Collars	
Ascend PD1 Pediatric Collar	SP30272-010
Ascend PD2 Pediatric Collar	SP30272-020
Ascend PD3 Pediatric Collar	SP30272-030
Ascend PD4 Pediatric Collar	SP30272-040
Ascend PD5 Pediatric Collar	SP30272-050
Back Panels	
Back Panel, Ascend Collar, Small	SP30372-020
Back Panel, Ascend Collar, Standard	SP30372-030
Back Panel, Ascend Collar, Large	SP30372-040
Front Panels	
Front Panel, Ascend Collar, Short	SP30472-020
Front Panel, Ascend Collar, Regular	SP30472-030
Front Panel, Ascend Collar, Tall	SP30472-040
Front Panel, Ascend Collar, X-Tall	SP30472-050
Pad Kits	
Pad Kit, Ascend Collar, Small & STND	SP30010-030
Pad Kit, Ascend Collar, Large	SP30010-040
Pad Kit, Ascend Collar, PD1/PD2	SP30020-030
Pad Kit, Ascend Collar, PD3/PD4/PD5	SP30020-040

